

**Innovatively  
designed** from  
foot to handgrip


**Crafted with care**  
from high-quality  
materials

**Competitive  
pricing** that  
improves your  
bottom line



# Big Foot Cane™

WORLD'S BEST CANE for **STABILITY**



## **BIG FOOT™ = STABILITY**

It all starts with our patented Big Foot. It's 4–5 times bigger than a standard cane tip, and articulates for exceptional stability. Front and rear shock-absorbing dampeners decrease stress from walking, and the foot's "grippy" tread offers superior traction.




## **CROSSBOW TECH™**

Patent-pending Crossbow Tech creates ultimate balance and control. The calculated alignment of the Big Foot, ergonomic handgrip, and 12° Crossbow bend creates a neutral balance point for the user. Result? Greater stability.



## **COMPRESSIVE HANDGRIP**

The comfortable, forward-facing handgrip reduces pressure on wrist and hands. Ideal for patients with arthritis or hand injuries.



## **ABOUT US**

Our innovative mobility devices keep people moving with dignity and vigor — and help them feel more in control. Invented by a world-class athlete who spent 1000+ days on traditional crutches, canes, and walkers.

## **GET IN TOUCH**

To learn more about our complete line of innovative mobility devices, contact [jmiller@iconmedical.biz](mailto:jmiller@iconmedical.biz)



Big Foot Cane



Hurrycane



Posture Cane



Quad Cane

Adjustable	✓	✓	✓	✗
Lightweight	✓	✓	✓	✗
Portable	✓	✓	✓	✗
Affordable	✓	✓	✓	✗
Easy-to-hold compressive handgrip	✓	✗	✗	✗
Sit-stand assist	✓	✗	✗	✗
Rattle-free use	✓	✗	✗	✗
All-terrain foot	✓	✗	✗	✗
Crossbow Tech	✓	✗	✗	✗
Approved by medical professionals	✓	✗	✗	✗


extends to fit 5'2" – 6'2"

SINGLE CANE


max. capacity 300 lbs.

PAIR OF CANES

max. capacity 450 lbs.



Ella  
BLACK



Oliver  
GREEN

“Most innovative cane I’ve seen in my career.”

— Dr. Howard M.  
Sports Medicine



**PRODUCT SPECIFICATIONS**  
Model No. V10ML — Big Foot Cane

HEIGHT	TUBE	HANDGRIP (W x D)	WEIGHT
32.3 in.	.75 in	1.75 x 4.75 in	1.3 lbs
82 cm	1.9 cm	4.44 x 12.06 cm	0.60 kg
BIG FOOT: 2.25 x 6 in / 5.71 x 15.24 cm			
COLORS: Ella Black (V10MLXAXX) / Oliver Green (V10MLMAXX)			

Products can be ordered with or without retail packaging.

[jmiller@iconmedical.biz](mailto:jmiller@iconmedical.biz)

**KEY PARTNERS**

We are partnering with all major retail, DME, and omni-channel marketers to develop a global footprint. We invite you to join us in providing mobility aids that can help people around the world.

**GET IN TOUCH**

To learn more about our complete line of innovative mobility devices, contact [jmiller@iconmedical.biz](mailto:jmiller@iconmedical.biz)