





User is 6'3" tall



Compact size — plus helpful aid to get up & down

Adjusts from cane height to trekking pole height



Ella BLACK



Oliver GREEN

Buy in pairs for use as trekking poles.

extends 33.8 to 56 inches  
fits users 4'6"–6'6"

SINGLE PRIMESTICK  
max. capacity 150 lbs

PAIR OF PRIMESTICKS  
max. capacity 250 lbs

“Multi-purpose device for active users, and a good-looking cane alternative. I highly recommend it!”

— Dr. Charles Mehlman



DYNAMO®

PRODUCT SPECIFICATIONS  
Model No. S1DML — PrimeStick

HEIGHT	TUBE	HANDGRIP	WEIGHT
37.3 in	.66 in	5 in (height)	1 lb
94.7 cm	1.67 cm	12.7 cm	.46 kg
CYCLONE FOOT: 2.25 in (diameter) / 15.24 cm (diameter)			
COLORS: Ella Black (S1DMLXAXX) / Oliver Green (S1DMLMAXX)			

Each PrimeStick comes in a handy mesh bag.

[jmiller@iconmedical.biz](mailto:jmiller@iconmedical.biz)

## GET IN TOUCH

To learn more about our complete line, contact [customerservice@medtecllc.com](mailto:customerservice@medtecllc.com)

## KEY PARTNERS

We are partnering with all major retail, DME, and omni-channel marketers to develop a global footprint. We invite you to join us in providing mobility aids that can help people around the world.